GO-ON® Product Information GB

Description GO-ON® is a sterile non-pyrogenic solution of Sodium Hyaluronate.

Composition A syringe of GO-ON® contains 2.5 ml of a 1% Sodium Hyaluronate solution as well as sodium chloride, sodium monohydrogen phosphate, sodium dihydrogen phosphate and water for injection.

Characteristics

Sodium Hyaluronate is a natural high molecular weight polysaccharide composed of a linear chain of disaccharide units made by sodium gluconate and N-acetylglucosamine. Sodium Hyaluronate is widely distributed throughout the tissues of the human body and is present in high concentrations in the synovial fluid, of which it is the major component. Sodium Hyaluronate plays a significant role in modulating the physical and mechanical interaction between adjacent tissues functioning as a tissue lubricant; it acts also as a viscoelastic support maintaining separation between tissues. Different Sodium Hyaluronate preparations may have different molecular weights but they have the same chemical structure. GO-ON® is a 1% Sodium Hyaluronate solution obtained from Streptococcus Equi by fermentation processes and subsequent purification.

Indications

GO-ON® is indicated as a viscoelastic supplementation for synovial fluid in knee and shoulder joints. GO-ON® is also indicated as viscoelastic supplementation in other synovial joints. The actions of the product are lubrication and mechanical support and it is suited for treatment of the symptoms of osteoarthritis.

Contraindications

GO-ON® must not be used in patients with a history of hypersensitivity to any of the ingredients of the product and in the case of inflammatory joint diseases such as rheumatoid arthritis or Bechterew disease.

Warnings and Precautions

Long-term viscoelastic supplementation with GO-ON® in synovial joints other than knee and shoulder joints has not been used only in cases appropriate for joint replacement.

This product should only be used by qualified persons who are familiar with its use.

The general precautions for intra-articular injections should be observed and, particularly, care should be taken to ensure that the site of injection is aseptic.

Preparation must be discarded in the event of damage to the ready-prepared syringe or sterile container. The product must not be used after the expiry date indicated on the box. The product must be kept out of the reach of children.

Interactions with Other Agents

There is a known incompatibility between sodium hyaluronate and quaternary ammonium salts such as benzalkonium chloride which can be used as skin disinfectants. GO-ON® must therefore never come into contact with such products, nor with medical or surgical equipment treated with these types of products, since hyaluronic acid can precipitate in their presence. To date there are no reports of incompatibilities (interactions) of GO-ON® with other solutions for intra-articular administration. The patient may benefit from the concurrent use of oral analgesic and anti-inflammatory agents during the first days of treatment.

Packaging

GO-ON® is available as a disposable sterile syringe containing a 2.5 ml solution of Sodium Hyaluronate (25 mg). Each syringe is contained in a sterile container. Packages contain 1, 3 or 5 syringes.

Shelf-life

Three years

Storage conditions

Store at a temperature lower than 25° C.

Symbols on the package

Anaphylactic reactions may occur. If this occurs, discontinue administration, closely monitor the patient and provide adequate treatment. The injection of GO-ON® may be followed by local side effects in the treated joint, e.g. pain, a sensation of heat, redness and/or swelling, effusion, infection. These adverse events can be reduced by the application of an ice-bag to the treated joint for 5 to 10 minutes.

Interactions with Other Agents

There is a known incompatibility between sodium hyaluronate and quaternary ammonium salts such as benzalkonium chloride which can be used as skin disinfectants. GO-ON® must therefore never come into contact with such products, nor with medical or surgical equipment treated with these types of products, since hyaluronic acid can precipitate in their presence. To date there are no reports of incompatibilities (interactions) of GO-ON® with other solutions for intra-articular administration. The patient may benefit from the concurrent use of oral analgesic and anti-inflammatory agents during the first days of treatment.

Packaging

GO-ON® is available as a disposable sterile syringe containing a 2.5 ml solution of Sodium Hyaluronate (25 mg). Each syringe is contained in a sterile container. Packages contain 1, 3 or 5 syringes.

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ROTTAPHARM Ltd.,